

Phonak

Power Pack

Technical Data

Dimensions (L x W x H)	84.8 x 59 x 22 mm (3.3 x 2.3 x 0.9 ")
Weight	37 g (1.3 oz)
Color	Champaign
Operating conditions	5° to +40° Celsius (41° to 104°Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation	Up to 90% (non condensing)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 85% (non condensing)
Atmospheric pressure during use	800 hPa to 1500 hPa
Included Parts	micro-USB Cable (for charging only)
External Interfaces	micro-USB port
Battery	Li-Polymer
Charging time	3 hours



Power supply

Input	Type: DC, stabilized 5 V +/- 5% Max. current: 500 mA The battery cannot be over-charged. It may become warm during charging. Important notice: It is recommended to only use the original Phonak Chargers and Power Pack
-------	---

